

## PRODUCT DATA SHEET

PRODUCT:	<b>Polyclonal Antibody to Bcl-X<sub>L</sub>/s</b> [Anti-BCL2L PAb]
PRODUCT NUMBER:	B-1225
LOT NUMBER:	A9073
FORMULATION:	Liquid. Purified antibody in 0.05M sodium phosphate buffer containing 0.1% sodium azide and 0.2% gelatin.
CONCENTRATION:	0.1 mg/ml
APPLICATIONS:	<p><u>Immunohistochemistry:</u></p> <p>Frozen sections: dilution 5µg/ml; no pre-treatment needed; mouse and rat frozen intestinal tissue was fixed with Zamboni's fixative prior to staining.</p> <p>Paraffin sections: dilution 2.5µg/ml; pressure cooker antigen retrieval is recommended.</p> <p><u>Western Blot:</u> Dilution 2.5µg/ml; chemiluminescent detection; positive control tested was MCF-7 cells stimulated for 72 hours with 0.2µg/ml doxorubicin; normal rabbit IgG as negative control. <b><i>Human samples only!</i></b> Antibody <u>titration</u> is recommended for optimal results in individual <u>systems</u>.</p>
SOURCE/HOST:	From rabbit

# A.G. Scientific, Inc.

6450 Lusk Blvd. Suite E102  
SAN DIEGO, CA  
92121



---

**IMMUNOGEN:** Synthetic peptide corresponding to aa 201-216 of human Bcl-X.

**SPECIFICITY:** Recognizes human, mouse, rat and opossum Bcl-X<sub>L</sub> (strongly) and Bcl-X<sub>s</sub> (moderately).

**STORAGE & HANDLING:** Store at +4°C. Do not freeze.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: 