

CERTIFICATE OF ANALYSIS

PRODUCT: **Caspase-9 Substrate (Chromogenic)**
[Ac-LEHD-pNA; Ac-Leu-Glu-His-Asp-pNA
(pNA = p-Nitroaniline)]; C₂₉H₃₈N₈O₁₁, M.W. 674.7

PRODUCT NUMBER: C-1177

LOT NUMBER: C1296

APPEARANCE: White powder

SOLUBILITY: Soluble in DMSO. Can be further diluted with buffer, pH 7.5.

PURITY (HPLC): > 98.9%

PEPTIDE CONTENT: 83%

IDENTITY (MS): Identical to standard reference

**AMINO ACID
ANALYSIS:**

SPECIFICATIONS:

Asx: 0.80-1.20

Glx: 0.80-1.20

His: 0.80-1.20

Leu: 0.80-1.20

RESULTS:

Asx: 1.00

Glx: 0.96

His: 1.01

Leu: 1.02

PRODUCT DESCRIPTION: Chromogenic substrate for caspase-9 and cysteine proteases. Similar to Ac-LEHD-AMC (Prod.# C-1179) and Ac-LEHD-AFC but cleavage is monitored colorimetrically by absorbance at 405nm. λ_{max} of pNA is 400nm.

PROTOCOL:

HEPES-Buffer: 100 mM HEPES, pH 7.5; 20% (v/v) glycerol; 5 mM DTT, 0.5 mM EDTA.

Substrate: Prepare 1 mM stock solution in DMSO.

- Induce apoptosis and prepare cell lysate or use recombinant caspase.

- Prepare reaction solution:

A.G. Scientific, Inc.

6450 Lusk Blvd. Suite E102

SAN DIEGO, CA 92121



30 μ l of substrate stock solution + 240 μ l HEPES-Buffer + 30 μ l of cell lysate.

- Incubate for 1 hour at 37 °C.
- Measure with microplate reader at 400nm.
- Suggested controls: Reaction mixture...
...without substrate.
...with non-apoptotic cell lysate.
...with apoptotic cell lysate and caspase inhibitor.

STORAGE & HANDLING:

Keep cool and dry at -20°C .

CAUTION: For laboratory research & scientific manufacturing use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified:

A handwritten signature in blue ink, appearing to be 'J. Smith'.