



## PRODUCT DATA SHEET

<b>PRODUCT:</b>	<b>Polyclonal Antibody to Caspase-3 (active)</b>
<b>PRODUCT NUMBER:</b>	C-1337
<b>LOT NUMBER:</b>	A8865
<b>FORMULATION:</b>	Liquid. Affinity purified antibody in PBS containing 50% glycerol, 1% BSA and 0.02 sodium azide.
<b>CONCENTRATION:</b>	0.2 mg/ml
<b>SOURCE / HOST:</b>	Rabbit
<b>IMMUNOGEN:</b>	Synthetic peptide around the cleavage site of human Caspase-3
<b>SPECIFICITY:</b>	Recognizes human, mouse and rat active caspase-3 (p20 subunit). Does not cross react with procaspase-3 (~32kDa)
<b>APPLICATION:</b>	Immunohistochemistry (paraffin sections (10-20µg)) Western Blot (1µg/ml) Optimal conditions must be determined individually for each application.
<b>STORAGE &amp; HANDLING:</b>	Store as supplied at -20°C. Avoid freeze/thaw cycles. After opening, prepare aliquots and store at -80°C.  CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified:

A handwritten signature in black ink, appearing to be 'JK'.