



PRODUCT DATA SHEET

PRODUCT:	Polyclonal Antibody to Caspase-9 (active)
PRODUCT NUMBER:	C-1355
LOT NUMBER:	A8927
SOURCE / HOST:	Rabbit
IMMUNOGEN:	Synthetic peptide around the cleavage site of human caspase-9.
SPECIFICITY:	Recognizes human and rat active caspase-9. Detects a band of ~ 37 kDa by Western Blot. Does not cross-react with procaspase-9.
APPLICATION:	Immunocytochemistry: 10-20 µg/ml Immunoprecipitation: 10-20 µg/ml Western Blot: 0.5-1 µg/ml *Optimal conditions must be determined individually for each application.
PURITY DETAIL:	Epitope –affinity purified.
CONCENTRATION:	0.2 mg/ml

A.G. Scientific, Inc.

6450 Lusk Blvd. Suite E102
SAN DIEGO, CA
92121



FORMULATION:

Liquid. In PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide.

STORAGE & HANDLING:

Keep cool and at -20°C. Avoid freeze/thaw cycles.
After opening prepare aliquots and store at -80°C.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified:

A handwritten signature in black ink, appearing to be 'JK'.