
CERTIFICATE of ANALYSIS

PRODUCT:	Calmodulin Separopore
PRODUCT NUMBER:	C-2877
LOT NUMBER:	H1424
APPEARANCE:	Clear solution
LIGAND:	Calmodulin
LIGAND COUPLING METHOD :	CNBr activation
SOURCE OF LIGAND:	Pure Bovine brain Calmodulin
LIGAND DENSITY:	Approx. 1 mg calmodulin/mL drained resin
Matrix:	Agarose beads, 4%
EXCLUSION LIMIT(MR):	Approx. 2×10^7
PARTICLE SIZE (RANGE):	53-180 μm
pH STABILITY:	4-9
STORAGE & HANDLING:	Store desiccated at +4°C. CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: SK