



## CERTIFICATE OF ANALYSIS

**PRODUCT:** PI3-Kinase SH3 Domain Binding Peptide  
[Arg-Lys-Leu-Pro-Pro-Arg-Pro-Arg-Arg]  
M.W. 1174.7

**PRODUCT NUMBER:** P-1143

**LOT NUMBER:** A9026

**APPEARANCE:** Lyophilized solid

**SOLUBILITY:** Soluble in H<sub>2</sub>O

**PURITY:** >98% (by HPLC)

**PEPTIDE CONTENT:** 81.2%

**SUPPLIED AS:** Trifluoroacetate salt; 1mg net peptide/vial

**SEQUENCE DESCRIPTION:** Sequence is derived from directed random library using phosphoinositol 3-kinase (PI 3-K) SH3 domain as the ligand.

**STORAGE & HANDLING:** Keep cool and dry at -20°C.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified:

A handwritten signature in black ink, appearing to be 'JK'.