



CERTIFICATE OF ANALYSIS

PRODUCT: Polyclonal Antibody to PPAR δ (mouse)

PRODUCT NUMBER: P-1199

LOT NUMBER: A8457

SOURCE / HOST: From rabbit

APPLICATIONS: Western Blot: 1:500
Immunohistochemistry (Fr): 1:200

FORMULATION: 100 μ l of epitope affinity-purified rabbit IgG
(1 mg/ml) containing 1 mg/ml BSA and 0.05% sodium azide.

IMMUNOGEN: Synthetic peptide: M(1)E Q P Q E E T P E A R E E(14)C

SPECIFICITY: By Western Blot, this antibody detects an ~49kDa protein representing PPAR δ from NIH-3T3-L1 cell lysate.

STORAGE & HANDLING: Storage at -20 °C is recommended. Avoid freeze/thaw cycles.

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: 