

## CERTIFICATE of ANALYSIS

**PRODUCT:** **Protein Solubilizer X-100**  
 [Active Ingredient: Supplied as a highly purified, Sterile, 10% aqueous solution of Triton X-100]

**PRODUCT NUMBER:** P-1500

**LOT NUMBER:** J1082A/355897

**DESCRIPTION:** Clear colorless solution filled in 5 ml amber colored glass ampoules

TEST	SPECIFICATIONS:	RESULTS:
TRITON X-100 CONCENTRATION:	10.0 ±0.5%	9.98%
ALDEHYDE CONCENTRATION:	≤0.1 mM	0.08 mM
PEROXIDE CONCENTRATION:	≤0.1 mM	0.06 mM
CONDUCTIVITY:	≤5 μMho	2.5 μMho

**STERILITY:** Complies to USP standards

**STORAGE & HANDLING:** Store at +20°C. ***PROTECT FROM LIGHT!***

***FOR LABORATORY USE ONLY!***