

## Certificate of Analysis

**PRODUCT:** **Protein Solubilizer X-100**  
[Active Ingredient: Supplied as a highly purified, Sterile, 10% aqueous solution of **Triton X-100**]

**PRODUCT NUMBER:** P-1500

**LOT NUMBER:** L1337/356272

**DESCRIPTION:** Clear colorless solution filled in 5 ml amber colored glass ampoules

TEST	SPECIFICATIONS:	RESULTS:
TRITON X-100 CONCENTRATION:	10.0 ±0.5%	9.96%
ALDEHYDE CONCENTRATION:	≤0.1 mM	0.07 mM
PEROXIDE CONCENTRATION:	≤0.1 mM	0.06 mM
CONDUCTIVITY:	≤5 μMho	2.0 μMho

**STERILITY:** Complies to USP standards

**STORAGE & HANDLING:** Store at +20°C. **PROTECT FROM LIGHT!**

**FOR LABORATORY USE ONLY!**