
CERTIFICATE OF ANALYSIS

PRODUCT:	Recombinant Protein A/G, Proteina/g, protein ag, protein a/g, proteinag, protein-ag, protein-a/g, proteinag, Protein A/G, rProtein A/G, r-protein A/G, recombinant protein A/G, IgG purification, antibody purification, Protein A, Protein G, Protein A&G, M.W. 59.7 kDa.
PRODUCT NUMBER:	P-2415
LOT NUMBER:	E1079
APPEARANCE:	Lyophilized.
SOURCE:	Recombinant expressed in E. coli.
PURITY: (By SDS-PAGE and HPLC Analysis)	≥ 98%
ENDOTOXIN:	<0.1 EU/mg (Protein A/G)
THEORETICAL PI:	4.89
SOLUBILITY:	Clear solution in 5 mg/ml H ₂ O
STORAGE & HANDLING:	Store desiccated at -20°C. Avoid freeze/thaw cycles, Centrifuge the vial prior to opening!

CAUTION: For research use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: AD