

Certificate of Analysis

PRODUCT:	PBS 1X Solution CAS# [7647-14-5], Sodium Chloride CAS# [7447-40-7], Potassium Chloride CAS# [7558-79-4], Sodium Phosphate dibasic CAS# [7778-77-0], Potassium Phosphate monobasic
PRODUCT NUMBER:	P-2426
LOT NUMBER:	L1444
APPEARANCE:	Clear and Colorless solution
pH:	7.50
DESCRIPTION:	Phosphate Buffered Saline is a buffer solution commonly used in biological research. PBS has many uses because it is isotonic and non-toxic to most cells. These uses include substance dilution and cell container rinsing.
STORAGE & HANDLING:	Store at Room Temperature (+20°C). CAUTION: For laboratory research & scientific manufacturing use only. Not for human or drug use. The pharmacological and toxicological properties of this product have not been fully investigated. Use caution when handling. Do not use this compound if you are not fully trained or are unaware of the hazards involved.

Verified: KLS