

BIODETERGENTS

BIODETERGENTS AT AG SCIENTIFIC

Biological detergents are a class of molecules that can manipulate hydrophobic and hydrophilic interactions among molecules in a biological sample. Also referred to as biode detergents, these products are used for the study of membrane proteins due to their ability to lyse cells, solubilize membrane proteins and lipids, prevent non-specific binding, and control protein crystallization.

AG Scientific is proud to offer a wide variety of biological detergents that enhance and accelerate disease research, gene studies, and drug development. We utilize the most advanced practices to ensure you receive high-purity biode detergents for all your life science applications.

“AG Scientific has been outstanding in supplying us the highest quality materials at reasonable prices with dependability and integrity.”

— Kenny Fink, Operations Manager, MTI-Globalstem (Part of ThermoFisher)

ThermoFisher
SCIENTIFIC

The world leader in serving science

CRITICAL CHEMICAL COMPONENTS, SIMPLIFIED SOURCING

- Streamlined procurement process
- On time, every time
- Import/export expertise
- Custom packaging and labeling capabilities
- OEM/Bulk quantities available

Our team has the know-how to ensure your orders are processed efficiently and accurately. Let us handle the sourcing while you focus on your research.

SUPERIOR QUALITY, CONSISTENT RESULTS

To ensure accurate results are consistently produced, all of our products are subject to rigorous quality assurance procedures. We have maintained an ISO certified facility since 2005.

YOUR ALLY IN DISCOVERY

Each of our superior-quality antibiotics are produced with the advancement of scientific research in mind. Benefits of partnering with AG Scientific include:

- Single point of contact
- Personalized service
- Superior quality guaranteed
- Vast sourcing network
- Reduced overhead
- Checks & balances (ISO Certified)
- Reliability: our products have been cited >3,000 times!

AG Scientific, Inc. was founded in 1997 with the idea that excellent product sourcing, competitive pricing, and high-quality service with a human touch was possible in the life science industry. We have the connections you need to procure critical components for all life science applications.

POPULAR BIODETERGENTS

PRODUCT#	PRODUCT NAME
C-1019	CHAPS
C-1020	CHAPS 10% Solution
D-1304	n-Dodecyl-Beta-D-Maltopyranoside
D-2425	Digitonin
N-1470	Nuclease Away for DNase and RNase removal
N-2366	NP-40 (10%) in Water
O-1036	n-Octyl-Beta-D-Glucopyranoside
S-1020	Sulfobetaine-10
S-2066	Sodium Dodecyl Sulfate
T-2835	TWEEN® 80, 10% Solution

DETERGENTS FOR ALL APPLICATIONS

We supply biodegradable detergents for all life science applications including common bioprocessing methods, gel electrophoresis, drug manufacturing, immunostaining, and DNA extraction.

- Ionic detergents
- Non-ionic detergents
- Zwitterionic detergents
- Non-detergent sulfobetaines (NDSB)

We exist to be your sourcing ally in limitless life science discovery. Don't see the product you need? Let us know and we'll find it for you!



TRUSTED BY INDUSTRY-LEADING SCIENTISTS

Our products have been cited more than 3,000 times by industry-leading researchers and organizations all around the globe.



OUR PLEDGE TO SERVE YOU:

- We will listen to your requirements and be extremely responsive to your needs. Comprehensive customer service, communication, and accountability will take place at all phases of your project.
- We take great pride returning to a true level of customer responsiveness that has become rare in our industry.
- We will focus on serving you with personalized interactions: You'll get a real person on the phone or video conference to handle your needs promptly.
- We will work to ensure you are 100% comfortable in our partnership.
- We will employ the highest quality standards to ensure the utmost accuracy, including care in our sourcing and packaging of your products with on-time delivery.
- You will be consistently informed of the status of your project so that you can confidently and effectively communicate with your organization.
- We are committed ensuring our service to you is the best overall experience and earn you as a life-long client.



PROTEIN SOLUBILIZER AMPOULES

PRODUCT#	PRODUCT NAME
P-1500	Protein Solubilizer X-100, Sterile 10% Triton X-100
P-1504	Protein Solubilizer 35, Sterile 10% BRIJ-35 Protein
P-1505	Solubilizer 40, Sterile 10% NP-40

